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SUMMARY AND EXPLANATION

PHE reagent is reagent used for generation of colour in test PHE in diagnostic kits. Detailed description of work with reagent PHE and interpretation of result of test is part of particular kit manual.

PRINCIPLE

Substrate L-phenylalanine is oxidatively deaminated with help of aminooxidase and flavoprotein to phenyl pyruvic acid. PHE reagent acts as chelating agent in reaction mixture causing formation of green colour.

DIAGNOSTIC KIT CONTAINS

- PHE reagent 20 ml (Ref. 3007)
- information leaflet

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- All samples and inoculated products has to be considered as potentially infectious and obvious safety regulations of your country must be respected while handling
- Do not use overexpired product.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

Kit contains irritable iron (III) chloride.

R 36/R38 Irritating to eyes and skin.

S 26 In case of contact with eyes rinse immediately with plenty of water and seek medical advice.

S 45 In case of accident or if you feel unwell seek medical advice immediately (show the label when possible).

STORAGE

PHE reagent can be stored in temperatures between +2 and +8°C. Expiration date is showed on label of each package.

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