COLI

Detekčný prúžok pre mikrobiológiu / Diagnostic strip for microbiology / Detekční proužek pro mikrobiologii / Mikrobiológiai tesztcsík

IVD

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Dátum poslednej revízie: 10. 3. 2014

**SUMMARY AND EXPLANATION**

COLI is a detection strip used for screening identification of *Escherichia coli*. Detection strip shows colour change by creating indigo colour on the base of metabolism of examining samples.

**PRINCIPLE**

Breakdown of substrate 4-nitrophenyl-b-D glucoronide is attended by colour change and by creating of yellow colouring of detection zone, this indicates positive result of GLR test.

Breakdown of substrate L-tryptophan is attended by colour change and by creating of pink colouring of detection zone, this indicates positive result of IND test.

**DIAGNOSTIC KIT CONTAINS:**

- 50 test strips for 50 determinations
- information leaflet

**WARNINGS AND SAFETY PRECAUTIONS**

- For *in vitro* diagnostics use.
- For professional use only.
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious safety valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

**STORAGE**

At dark place at +2 to +8 °C. Expiration date is placed on the label of each package.

**REQUIRED MATERIALS**

IND reagent (Ref. 3002)  
Tube with 0,5 ml of saline solution

**RECOMMENDED PROCEDURE:**

1. Prepare thick bacterial suspension ( $\geq 3$  McF) into the tube with saline.
2. Put in diagnostic strip COLI.
3. Incubate in thermostat by  $35 \pm 2^\circ\text{C}$ .
4. Evaluate result of GLR test after 2 hours. Extend incubation for 4h in case of negative result.
5. Add 2 drops of IND reagent after evaluation of GLR test and evaluate result of IND test after 5 minutes.

**INTERPRETATION OF REACTION**

**Positive result of tests GLR (yellow) and IND (pink) indicate *Escherichia coli*.**  
For the need of own proof function, use following recommended bacterial strains:

Bacterial strain	CCM / ATCC	GLR	IND
<i>Escherichia coli</i>	4225 / 35218	yellow	pink

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

**TIP:** Use negative control – follow recommended procedure, but without bacterial culture.

Kit components doesn't contain hazardous substances.

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