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Dátum poslednej revízie: 28. 5. 2014

SUMMARY AND EXPLANATION

TRB is detection strip used for differentiation of bacteria belonging to genus *Moraxella* a *Neisseria*. Detection strip changes colour to yellow on base of metabolism of examined sample.

PRINCIPLE

Breakdown of substrate by enzymatic hydrolysis is attended by colour change and creation of yellow colouring, this indicates positive result of TRB test.

DIAGNOSTIC KIT CONTAINS:

- 50 test strips for 50 determinations
- information leaflet

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious safety valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

STORAGE

At dark place at +2 to +8 °C. Expiration date is placed on the label of each package.

REQUIRED MATERIALS

Tube with 0,5 ml of saline solution

RECOMMENDED PROCEDURE:

1. Prepare thick bacterial suspension (≥ 2 McF) from fresh 24 hour culture.
2. Put in diagnostic strip and dip zone into suspension.
3. Incubate in thermostat by temperature 35 ± 2°C for 2 - 4 hours.
4. Evaluate result of TRB test.

INTERPRETATION OF REACTION

Yellow colouring indicates positive result of TRB test.

Bacterial strain	TRB
<i>Moraxella catarrhalis</i>	+ (yellow)
<i>Neisseria sp.</i>	- (pink)

For the need of own proof function, use following recommended bacterial strains:

Control strain	CCM / ATCC	TRB
<i>Moraxella catarrhalis</i>	4391 / -	+ (yellow)
clear water	-	- (pink)

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

Kit components doesn't contain hazardous substances.

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