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Dátum poslednej revízie: 10.03.2014

SUMMARY AND EXPLANATION

DMACA reagent is reagent used for generation of colour in test IND. The reagent is also suitable for anaerobic microorganisms for its great sensitivity.

PRINCIPLE

By the action of the enzyme tryptophanase indole is created and detected with DMACA reagent.

DIAGNOSTIC KIT CONTAINS

- DMACA reagent 20 ml (Ref. 3009)
- Information leaflet

WARNINGS AND SECURITY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated reagent liquidate according to valid regulations. Packaging materials are recyclable.

R 68/20/21/22 Health harmful - possible risk of irreversible effects by inhalation, contact with skin or if swallowed.
S 45 In case of accident or if you feel unwell seek medical advice immediately (show the label when possible).

STORAGE

Store DMACA reagent in dark at temperatures between +2 and +8°C. Expiration date is showed on label of each package.

REQUIRED MATERIALS

Filtration paper

RECOMMENDED PROCEDURE

1. Use pure 24 hour bacterial culture. Smear 2-3 well separated colonies on filtration paper. Isolate colonies from suitable cultivation medium (for example: blood agar)
2. Add 1 drop of DMACA
3. Read the reaction visually after 0.5-2 minutes.

Bacterial strain	CCM / ATCC	POSITIVE	NEGATIVE
<i>Escherichia coli</i>	4225 / 35218	green-blue	without colour change
Distilled water	-	without colour change	without colour change

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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