

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Control Serum Anti-Salmonella Typhimurium
for Widal reaction**

Catalogue No. (Package size) TS 1624 (1 ml)

Lot No. **88 09 10**

Expiry date **2015-09**

Intended use The control serum reacts in the Widal test with the homologous
Salmonella Typhimurium Test Antigen (SIFIN) with a titre of
1:1600 (tube) and 1:3200 (MTP). Deviations of \pm one titre step are
permitted.

Specificity The serum control (antigen addition is substituted by physiological
saline solution) must not show any agglutination.

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin, 01.10.10


Head of Quality Control