

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Control Serum Anti-Salmonella Paratyphi A</b> for the Widal Reaction
Catalogue no.	TS 1626
Package size	1 ml
<b>Lot no.</b>	<b>7780917</b>
<b>Expiry date</b>	<b>2022-09</b>
Intended use	The control serum agglutinates in the Widal test with the homologous Salmonella Paratyphi A-O Test Antigen (sifin) with a titre of 1:3200 (tube) and 1:6400 (MTP) and with the homologous Salmonella Paratyphi A-OH Test Antigen (sifin) with a titre of 1:3200 (tube) and 1:12800 (MTP). Deviations of $\pm$ one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by phosphate buffered saline solution) must not show any agglutination.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
place, date

22.09.2017

  
signature Head of Quality Control