

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321 13088 Berlin, Germany

T: +49 30 927 030-0 F: +49 30 927 030-30 info@sifin.de www.sifin.de

Product

Salmonella H:b Test Antigen

for the Widal Reaction

Catalogue no

TS 1632

Package size

10 ml

Lot no.

04 01 15

Expiry date

2016-01

Intended use

The Test antigen reacts in the Widal test with the homologous Salmonella antibody (Apti Hib) in the sorum to be tested. The specific bonding is indicated

antibody (Anti-H:b) in the serum to be tested. The specific bonding is indicated

by agglutination.

Specificity

There must be no reaction with physiological saline solution.

Preservatives

Formalin 0.5 %

Storage

2 °C...8 °C

Berlin, 26.03.41

place, date

signature Head of Quality Control