

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Salmonella H:c Test Antigen for the Widal Reaction
Catalogue no.	TS 1633
Package size	10 ml
Lot no.	64 03 15
Expiry date	2016-03
Intended use	The Test antigen reacts in the Widal test with the homologous Salmonella antibody (Anti-H:c) in the serum to be tested. The specific bonding is indicated by agglutination.
Specificity	There must be no reaction with physiological saline solution.
Preservatives	Formalin 0.5 %
Storage	2 °C...8 °C

Berlin, 16.04.15
 place, date


 signature Head of Quality Control