



## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Control Serum Anti-Salmonella H:d for the Widal Reaction</b>
Catalogue no.	TS 1644
Package size	1 ml
Lot no.	09 01 14L
Expiry date	2017-01
Intended use	The control serum agglutinates in the Widal test with the homologous Salmonella H:d Test Antigen (sifin) with a titre of 1:3200 (tube). Deviations of $\pm$ one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by phosphate buffered saline solution) must not show any agglutination.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 14.07.14  
place, date

  
signature Head of Quality Control