

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Control Serum Anti-Salmonella H d for Widal reaction</b>
Catalogue No. (Package size)	TS 1644 (1 ml)
Lot No.	<b>105 10 09</b>
Expiry date	<b>2014-09</b>
Intended use	The control serum reacts in the Widal test (tube) with the homologous Salmonella Test Antigen H d (SIFIN) with a titre of 1:3200. Deviations of $\pm$ one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by physiological saline solution) must not show any agglutination.
Preservatives	Sodlum azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 03.11.09

  
Head of Quality Control