

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Yersinia enterocolitica-O9 Test Antigen
for the Micro Widal Reaction**

Catalogue No.
(Package size) TS 1723 (3 ml)

Lot No. **54 02 13**

Expiry date **2018-06**

Intended use In the micro Widal test, the test antigen reacts with the homologous *Yersinia enterocolitica* O9 antibodies in the serum to be tested. The specific bonding is indicated by agglutination.

Specificity There must be no reaction with heterologous antibodies anti-*Yersinia* O3, O5 and O6 or with physiological saline solution.

Preservatives Thiomersal 0.2 mg/ml

Storage 2...8 °C

Berlin,

20.02.14



Head of Quality Control