



Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Yersinia enterocolitica-O 9 Test Antigen for the Widal Reaction
Catalogue no.	TS 1723
Package size	3 ml
Lot no.	71 06 13L
Expiry date	2018-06
Intended use	in the Widal test (microtitre plate) the dyed test antigen reacts with the specific <i>Yersinia enterocolitica</i> O 9 antibodies in the serum to be tested. The specific bonding is indicated by agglutination.
Specificity	There must be no reaction with nonspecific antibodies Anti- <i>Yersinia enterocolitica</i> O 3, O 5 and O 6 or with physiological saline solution.
Preservatives	Thiomersal 0.2 mg/ml
Storage	2 °C...8 °C

Berlin, 26.05.14
place, date


signature Head of Quality Control