

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh


Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Yersinia enterocolitica-O 9 Test Antigen for the Widal Reaction
Catalogue no.	TS 1723
Package size	3 ml
Lot no.	8320918
Expiry date	2023-09
Intended use	In the Widal test (microtitre plate) the dyed test antigen reacts with the specific <i>Yersinia enterocolitica</i> O 9 antibodies in the serum to be tested. The specific bonding is indicated by agglutination.
Specificity	There must be no reaction with nonspecific antibodies Anti- <i>Yersinia enterocolitica</i> O 3 and O 5 or with physiological saline solution.
Preservatives	Thiomersal 0.2 mg/ml
Storage	2 °C...8 °C

Berlin,
place, date

27.02.2019


signature Head of Quality Control