

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Control Serum Anti-Yersinia enterocolitica O 3 for the Widal Reaction
Catalogue No. (Package size)	TS 1733 (1 ml)
Lot No.	69 04 12
Expiry date	2017-04
Intended use	The control serum agglutinates in the tube test with the uncoloured test antigen Yersinia enterocolitica O 3 (SIFIN) with a titre of 1:640; on the MTP it reacts with the coloured homologous test antigen (SIIFN) with a titre of 1:1600. Deviations of \pm one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by physiological saline solution) must not show any agglutination.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

16.05.12


Head of Quality Control