

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Control Serum Anti-Yersinia enterocolitica O 3 for the Widal Reaction
Catalogue no.	TS 1733
Package size	1 ml
Lot no.	68 04 16
Expiry date	2021-04
Intended use	The control serum agglutinates on the microtitre plate (MTP) with the coloured test antigen Yersinia enterocolitica O 3 (sifin) with a titre of 1:1600. Deviations of \pm one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by phosphate buffered saline or physiological saline solution) must not show any agglutination.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 18.05.2016

 place, date

J. Ludwig

 signature Head of Quality Control