

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Control Serum Anti-Yersinia enterocolitica O 5 for the Widal Reaction
Catalogue no.	TS 1734
Package size	1 ml
Lot no.	120 05 15
Expiry date	2020-05
Intended use	The control serum agglutinates on the microtitre plate (MTP) with the coloured test antigen Yersinia enterocolitica O 5 (sifin) with a titre of 1:1600. Deviations of \pm one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by phosphate buffered saline or physiological saline solution) must not show any agglutination.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,
 place, date

24.06.15

signature Head of Quality Control

