


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test Serum Anti-Shigella II (<i>S. dysenteriae</i> type 1-10)
Catalogue No. (Package size)	TS 1821 (1 ml)
Lot No.	146 11 11
Expiry date	2016-11
Intended use	Test Serum Anti-Shigella II records the serotypes of <i>Shigella dysenteriae</i> (type 1 to 10) in slide agglutination.
Specificity	No agglutination with the serotypes of <i>Shigella flexneri</i> , <i>Shigella boydii</i> (type 1 to 15) and <i>Shigella sonnei</i> S (Phase I)- and F-form (Phase II).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 19.12.12


Head of Quality Control