

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485 sifin diagnostics gmbh

Berliner Allee 317-321 13088 Berlin, Germany

T: +49 30 927 030-0 F: +49 30 927 030-30 info@sifin.de www.sifin.de

Product

Anti-Coli O 114:K 90

Catalogue no.

TS 2208

Package size

1 ml

Lot no.

8140918

Expiry date

2023-09

Intended use

Anti-Coli O 114:K 90 records all E. coli strains that contain the antigen O 114:K 90 in slide agglutination. In the Widal reaction the serum reaches a titre of 1:800 (tube and MTP) with the homologous E. coli strains; deviations of

±1 titre step are acceptable.

Specificity

Since unspecific co-agglutinations with other E. coli strains often occur in the slide agglutination test, a positive result in the slide test must always be confirmed by examination of the boiled E. coli isolate in the Widal test.

Preservatives

Sodium azide 0.9 mg/ ml

Storage

2 °C...8 °C

22,10.2018

place, date