

## PRS-4/12, platform

### DESCRIPTION

Can accommodate 4 tubes, 20-30 mm diameter (up to 50ml tubes) and 12 tubes, 10-16 mm diameter (1.5 ml - 15 ml tubes).

### SPECIFICATIONS

Tube diameter	20-30 and 10-16 mm
Tube volume	up to 50 and 1.5-15 ml
Capacity	4 and 12
Platform type	optional



### CAT. NUMBER

BS-010117-AK

Platform PRS-4/12



**Bio RS-24**  
Including PRS-22 platform  
Mini-rotator

Mini-rotator **Bio RS-24** provides vertical rotation of the platform. The rotator is an ideal instrument for preventing blood coagulation in tubes and for fulfilment of procedures of biological components extraction. ...

[read more](#)