

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

|                    |  |
|--------------------|--|
| <b>Product</b>     | <b>Anti-Coli O 25:K 11</b>   |
| Catalogue no.      | TS2201   |
| Package size       | 1 ml   |
| <b>Lot no.</b>     | <b>0990123</b>   |
| <b>Expiry date</b> | <b>2025-01</b>   |
| Intended use       | Anti-Coli O 25:K 11 records all <i>E. coli</i> strains that contain the antigen O 25:K 11 in slide agglutination. In the Widal reaction the serum reaches a titre of 1:1600 (tube) and 1:1600 (MTP) with the homologous <i>E. coli</i> strains; deviations of $\pm 1$ titre step are acceptable. |
| Specificity        | Since unspecific co-agglutinations with other <i>E. coli</i> strains often occur in the slide agglutination test, a positive result in the slide test must always be confirmed by examination of the boiled <i>E. coli</i> isolate in the Widal test.  |
| Preservatives      | Sodium azide 0.9 mg/ ml  |
| Storage            | 2 °C...8 °C  |

**This document has been established electronically and is valid without signature.  
The batch release was done in accordance with SOP QMV 7.**