

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella H:h</b>
Catalogue no.	TR1409
Package size	1 ml
<b>Lot no.</b>	<b>6750922</b>
<b>Expiry date</b>	<b>2027-08</b>
Intended use	Anti-Salmonella H:h records all <i>Salmonella</i> strains, as long as they contain antigen H:h (in the antigen combination H:e,h) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 24.10.2022  
place, date

  
signature Head of Quality Control