

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella H:z</b>
Catalogue no.	TR1424
Package size	1 ml
<b>Lot no.</b>	<b>0600123</b>
<b>Expiry date</b>	<b>2026-09</b>
Intended use	Anti-Salmonella H:z records all <i>Salmonella</i> strains in Subspecies I, II,III, as long as they contain antigen H:z in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.  
The batch release was done in accordance with SOP QMV 7.**