


### CERTIFICATE OF ANALYSIS

<b>Product Name:</b> Prolex™ <i>E. coli</i> O157 Latex Test Reagent Kit		<b>Lot Number:</b> C21745
<b>Product Code:</b> PL.070B		<b>Expiration Date:</b> 2025-03-31
<b>Product Code</b>	<b>Quantity</b>	<b>Components</b>
PL.072B	1 x 3.1 ml	<i>E. coli</i> O157 Latex Reagent
PL.074B	1 x 1.5 ml	<i>E. coli</i> O157 Positive Control
PL.077B	1 x 1.5 ml	<i>E. coli</i> O157 Latex Reagent Negative Control
	75	Mixing Sticks
PL.092-20	20 cards	Test Cards
<b>Additional Information</b>		
Storage Conditions:		2-8°C
Instructions for Use:		Refer to pro-lab.com/resources for current version
<b>Quality Control Testing</b>		
<b>Criteria</b>		<b>Result</b>
Strong agglutination of Latex Reagent when tested against <i>E. coli</i> serotype O157 isolates.		PASS
No agglutination reaction of Latex Reagent when tested against non O157 serotypes <i>E. coli</i> , and <i>E. hermanii</i> isolates.		PASS
No agglutination reaction of Negative Control when tested against <i>E. coli</i> and <i>E. hermanii</i> isolates.		PASS
Strong agglutination of the <i>E. coli</i> O157 Positive Control when tested against Latex Reagent.		PASS
No agglutination of the <i>E. coli</i> O157 Positive Control when tested against Negative Control.		PASS
<b>Release</b>		
2023-07-06		
Date	Quality Control	

Canada ☎ 800 268 2341 Fax 905 731 0206  
 20 Mural Street, Unit #4, Richmond Hill, Ontario,  
 L4B 1K3  
 support@pro-lab.com

U.S.A. ☎ (512) 832 9145 Fax (512) 832 6424  
 21 Cypress Blvd., Suite 1155, Round Rock,  
 Texas 78665-1046  
 support@pro-lab.us

U.K. ☎ 0151 353 1613 Fax 0151 353 1614  
 3 Bassendale Road, Bromborough, Wirral,  
 Merseyside CH62 3QL  
 uksupport@pro-lab.co.uk