



[www.pro-lab.co.uk](http://www.pro-lab.co.uk)

### Certificate of analysis

<b>Product Name:</b> Grams Staining Kit (Dilute Carbol Fuchsin)	<b>Lot number:</b> K37753
<b>Product Code:</b> PL.8057/25	<b>Expiry date:</b> 2025 02 28

Product Code	Quantity	Product description
PL.7000/25	1 X 250ml	Crystal Violet
PL.7003/25	1 X 250ml	Grams Iodine
PL.7006/25	1 X 250ml	Grams Differentiator
PL.7015/25	1 X 250ml	Dilute Carbol Fuchsin

#### Additional Information

Description: Violet, brown, clear and magenta solutions.	Concentrate Lot numbers: K37503, K37686, K37667	Storage conditions: 15-25°C
Production date: 2024/03/12	Instructions for Use:	Precautions: refer to SDS
		Go to <a href="http://pro-lab.co.uk">pro-lab.co.uk</a> for current version.

#### Quality Control Testing

Challenged with test organisms:  
*Escherichia coli* NCTC 12241  
*Oxford Staphylococcus aureus* NCTC 12981

Test results:  
 Conforms to performance specifications when used in standard Grams stain method, giving crisp staining of Gram positive and Gram negative bacteria with a clean counterstained background.

Comments:  
 We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

#### Release

2024-03-13	
Date	Signed

T: 0151 353 1613  
 3 Bassendale Road, Bromborough,  
 Wirral, CH62 3QL, UK.  
[uksupport@pro-lab.com](mailto:uksupport@pro-lab.com)