

Quality control certificate



Biochemical Wedge Agar Base

N° Cat: 1429

N° Batch: 303045

Re-test Date: 3/2026

Formula in g/l

Bromthymol blue	0,06	Glucose	6
Bacteriological agar	13	Beef extract	4
Peptone	12	Sodium chloride	5
Urea	6		

Typical formula g/L * Adjusted and/or supplemented as required to meet performance criteria.

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Batch pH
w/o rests	Fine powder	Beige	Green	6,6-7,0	6,6

Microbiological test

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.