

## **CERTIFICATE OF ANALYSIS**

Product Name:	Staph Latex Kit			Lot Number:	C22464
Product Code:	PL.080B			Expiration Date:	2025-08-31
Product Code	Quantity	Compone	Components		
PL.083B	2 x 2.5 ml	Staph Late	Staph Latex Reagent		
PL.085B	1 x 2.5 ml	Negative C	Negative Control Latex Reagent		
	200	Mixing Stic			
PL.092-20	20 cards	Test Cards	ards		
Additional Inform					
Storage Conditions: 2-8°C					
Instructions for Use: Refer to pro-lab.com/resources for current version					
Quality Control T	esting				
Criteria					Result
Strong agglutination of Staph Latex Reagent when tested against Staphylococcus aureus.					PASS
No agglutination of Negative Control when tested against test isolates.					PASS
Release					
2024-03-07			l	S.	
Date				Quality Control	

Canada 🕾 800 268 2341 Fax 905 731 0206 20 Mural Street, Unit #4, Richmond Hill, Ontario, L4B 1K3 support@pro-lab.com U.S.A. <sup>@</sup> (512) 832 9145 Fax (512) 832 6424 21 Cypress Blvd., Suite 1155, Round Rock, Texas 78665-1046 support@pro-lab.us U.K. 🕾 0151 353 1613 Fax 0151 353 1614 3 Bassendale Road, Bromborough, Wirral, Merseyside CH62 3QL uksupport@pro-lab.co.uk