

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella OMA O-Group-Pool: A, B, D, E, L
Catalogue no.	TR1151
Package size	1 ml
Lot no.	5900723
Expiry date	2028-06
Intended use	Anti-Salmonella OMA records all <i>Salmonella</i> strains belonging to the O-Groups A, B, D, E, L in slide agglutination.
Specificity	No agglutination with <i>Salmonella</i> strains which are recorded by Anti-Salmonella OMB (O-Group-Pool: C, F, G, H), Anti-Salmonella OMC (O-Group-Pool: I, J, K, M, N, O, P), Anti-Salmonella OMD (O-Group-Pool: Q, R, S, T, U, V, W), Anti-Salmonella OME (O-Group-Pool: X, Y, Z, 51, 52, 53), Anti-Salmonella OMF (O-Group-Pool: 54-59) and Anti-Salmonella OMG (O-Group-Pool: 60-63, 65-67).
Preservatives	Sodium azide 0,9 mg/ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**