

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella O:10
Catalogue no.	TR1308
Package size	1 ml
Lot no.	0830123
Expiry date	2027-12
Intended use	Anti-Salmonella O:10 records all <i>Salmonella</i> strains of groups E ₁ and E ₄ , as long as they contain antigen O:10 in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**