

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Coli I
Catalogue no.	TS2111
Package size	1 ml
Lot no.	5780723
Expiry date	2028-06
Intended use	Anti-Coli I records all <i>E. coli</i> strains, as long as they contain the antigens O 26:(K 60), O 44:(K 74), O 114:(K 90), O 125:(K 70), O 142:(K 86) and O 158:(K-).
Specificity	No agglutination with <i>E. coli</i> strains which are recorded by Anti-Coli II and Anti-Coli III.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**