


### CERTIFICATE OF ANALYSIS

|  |  |  |
|--|--|--|
| <b>Product Name: Staph Latex Kit</b>   |  | <b>Lot Number: C22831</b>                          |
| <b>Product Code: PL.080B</b>   |  | <b>Expiration Date: 2025-10-31</b>                 |
| <b>Product Code</b>  | <b>Quantity</b>  | <b>Components</b>                                  |
| PL.083B  | 2 x 2.5 ml   | Staph Latex Reagent                                |
| PL.085B  | 1 x 2.5 ml   | Negative Control Latex Reagent                     |
|  | 200  | Mixing Sticks                                      |
| PL.092-20  | 20 cards   | Test Cards   |
| <b>Additional Information</b>  |  |  |
| Storage Conditions:  |  | 2-8°C  |
| Instructions for Use:  |  | Refer to pro-lab.com/resources for current version |
| <b>Quality Control Testing</b>   |  |  |
| <b>Criteria</b>  |  | <b>Result</b>                                      |
| Strong agglutination of Staph Latex Reagent when tested against <i>Staphylococcus aureus</i> . |  | PASS   |
| No agglutination of Negative Control when tested against test isolates.                        |  | PASS   |
| <b>Release</b>   |  |  |
| 2024-07-04   |  |  |
| Date   | Quality Control  |  |

Canada ☎ 800 268 2341 Fax 905 731 0206  
 20 Mural Street, Unit #4, Richmond Hill, Ontario,  
 L4B 1K3  
 support@pro-lab.com

U.S.A. ☎ (512) 832 9145 Fax (512) 832 6424  
 21 Cypress Blvd., Suite 1155, Round Rock,  
 Texas 78665-1046  
 support@pro-lab.us

U.K. ☎ 0151 353 1613 Fax 0151 353 1614  
 3 Bassendale Road, Bromborough, Wirral,  
 Merseyside CH62 3QL  
 uksupport@pro-lab.co.uk