

## QUALITY CONTROL CERTIFICATE

### **CLOTRIMAZOLE 10 µg Neo-Sensitabs**

REF: 81212

LOT: 24-1027

Date of analysis: 2024-10

Date of expiry: 2029-10

Appearance: Light yellow, plane tablet, 9 mm in diameter,  
on both sides printed with a black identification  
code.  
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.  
Complies.

Performance:

Test strains:  
Cand. tropicalis ATCC 750  
C. albicans ATCC 10231

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets  
vary less than  $\pm 1$  mm from mea  
Complies.

Microbiological quality: Complies.

Rosco Diagnostica ApS



Claus Bo Svendsen MD PhD MBA  
Chief Medical Officer, Head of Quality