

Specification

A sterile selective supplement for the isolation of *Listeria* species.

Presentation

10 Freeze dried vials
Vial
with: 9 ± 0.1 g

Packaging Details

23x60 mm glass vials, tag labelled, White plastic cap - 10 vials per box.

Shelf Life

49 months

Storage

2-25 °C

Composition

Compositon (g/vial)

Sodium Nalidixate.....0.0100
Acriflavine.....0.0125
Ferric Ammonium Citrate.....0.2500

NOTE : Each vial is sufficient to supplement 500 ml of medium Base: Fraser Borth Base.

Reconstitute the original freeze-dried vial
by adding
Sterile Distilled Water..... 6 ml

Description /Technique

Description:

This supplement is added in Fraser broth base in order to obtain a secondary enrichment complete medium.

This medium is a modification of the UVM broth. It gives better results in the detection rate of *Listeria monocytogenes* in meat products and has the added advantage of only taking 3-4 days.

Technique:

Collect, dilute and prepare samples and volumes as required according to specifications, directives, official standard regulations and/or expected results.

Reconstitute the vial with 6 ml of the sterile diluent in aseptic conditions and add it to 500 ml of sterilized Broth base cooled to 50°C.

Do not overheat once supplemented.

Pour the complete medium into tubes and inoculate.

Incubate the tubes in aerobic atmosphere at 37 ± 1°C for 24±2h.

Incubation times longer than those mentioned above or different incubation temperatures may be required depending on the sample or the specifications.

After incubation, the isolation is carried out on *Listeria* agar according to Ottaviani & Agosti and a second selective agar for *Listeria*, eg Oxford, Palcam, or any other selective agar.

In these media, the colonies that present blackening due to the hydrolysis of esculin are presumptively typical strains of *Listeria*.

Quality control

Physical/Chemical control

Color : Yellowish-brown

pH: at 25°C

Microbiological control

Add 1 vial to 500 ml of medium base. DO NOT HEAT once supplemented.

Inoculate 30-300 CFU (productivity) 1.000-10.000 CFU (selectivity)

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

Aerobiosis. Incubation at 35°C ± 2 °C, reading at 24-48 hours

Microorganism

L. monocytogenes ATCC® 13932, WDCM 00021

Escherichia coli ATCC® 25922, WDCM 00013

Listeria monocytogenes ATCC® 35152

Sterility Control

Add 5mL of the sample to 100 mL of TSB and to 100 mL Thioglycollate.

Incubation 48 h at 30-35 °C and 48 h at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

Growth

Good. Black medium. Positive esculine

Inhibited. Confirm in TSA at 37°C±1 reading 24 ± 3h

Good. Black medium. Positive esculine

Bibliography

- ATLAS, R.M. (1993) Handbook of Microbiological Media. CRC Press. Boca Raton. Florida.
- FRASER, J.A. & W.H. SPERBER (1988) Rapid detection of *Listeria* spp. In food and environmental samples by esculin hydrolysis. J. Food Prot. 51:762-765.
- ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- ISO 11290-1:2017 Standard. Microbiology of the food chain. Horizontal method for the detection and enumeration of *Listeria monocytogenes* and for *Listeria* spp.- Part 1: Detection Method
- ISO 11290-2:2017 Standard. Microbiology of the food chain. Horizontal method for the detection and enumeration of *Listeria monocytogenes* and for *Listeria* spp.- Part 2: Enumeration Method.
- McCLAIN, D. & W.H. LEE (1988) Development of a USDA-FSIS method for isolation of *Listeria monocytogenes* from raw meat and poultry. J.AOAC 71:660-664.
- VANDERZANT, C & D.F. SPLITTSTOESSER (1992) Compendium of methods for the microbiological examination of foods. APHA. Washington. DC.