

Specification

Selective supplement for the isolation of *Burkholderia cepacia*.

Presentation

10 Freeze-dried vials
Vial
with: 6 ± 0.1 ml

Packaging Details

22±0.25 x 55±0.5 mm glass vials, tag labelled, White plastic cap - 10 vials per box.

Shelf Life

49 months

Storage

2-25 °C

Composition

Composition (g/vial):

Polymyxin B Sulphate..... 75.000 IU
Gentamicin 0.0025
Ticarcillin..... 0.0500

NOTE: Each vial is sufficient to supplement

500ml of medium Base for *Burkholderia cepacia* spp.

Reconstitute the original
freeze-dried vial
by adding 1 vial with sterile
distilled water..... 5 ml

Description /Technique

Description:

Burkholderia Cepacia Agar Base (Cat. 1347) is a selective medium especially formulated for the isolation of *Burkholderia cepacia* (*Pseudomonas cepacia*), from clinical and non-clinical specimens. *Burkholderia cepacia* is a Gram-negative, oxidase positive, mobile and aerobic bacillus. It is normally found in water deposits and damp environments. This bacillus is an important opportunist pathogen and causes pulmonary infections in Cystic Fibrosis patients.

The organism may be present in small numbers in many non sterile products used in hospitals. It has been isolated from a number of water sources and can grow in distilled water with a nitrogen source because of its capacity to fix CO₂ from air. Suction catheters rinsed in a solution of acetic acid have reduced the transmission of *Burkholderia cepacia* and other *Pseudomonas*.

The medium contains peptone, which provides nitrogen, vitamins, minerals and amino acids essential for growth. Selective agents are added to improve *B. cepacia* recovery through the inhibition of common contaminants. Crystal violet inhibits Gram-positive cocci, especially enterococci and staphylococci. Bile salts inhibit most Gram-positive cocci except for enterococci, and ticarcillin and polymyxin B inhibit Gram-negative bacilli. Phenol red facilitates detection of *B. cepacia*. Alkaline end products from the metabolism of pyruvate raise the pH of the medium, causing the color of the indicator to change from light orange to pink, or pink-red, in the growth area. In areas of heavy *B. cepacia* growth, the pink color intensifies. Magnesium sulphate, ammonium sulphate and ferrous sulfate provide sources of sulfates and metallic ions. Phosphate salts act as a buffer system. Bacteriological agar is the solidifying agent.

Colonies of *B. cepacia* are 1-2 mm in diameter and turn the medium to pink. Low numbers of colonies may not produce a color change of the medium. Occasional growth of some strains of *Candida* species, *Stenotrophomonas maltophilia*, *Pseudomonas aeruginosa* and other *Pseudomonas* species may occur.

Technique:

Aseptically reconstitute 1 vial with 5 ml of sterile distilled water. Mix gently until complete dissolution and aseptically add to 500 ml of *Burkholderia Cepacia* Agar Base (Cat. 1347), autoclaved and cooled to 50 °C. Mix well and distribute into sterile containers.

Instructions for use:

- Take a routine respiratory sample from the patient (e.g. sputa, deep pharyngeal swabs or bronchial washings).
- Dilute the sample if necessary.
- Streak onto *Burkholderia cepacia* Agar Base and incubate at 37 °C for 48-72 hours.
- Examine for green colonies and the medium turning to bright pink.
- Re-incubate for a further 24 hours if necessary.

Quality control

Physical/Chemical control

Color : White-Gray

pH: at 25°C

Microbiological control

Reconstitute 1 vial as indicated in COMPOSITION; shake and dissolve completely

Add 1 vial to 500 ml of medium base. DO NOT HEAT once supplemented.

Aerobiosis. Incubation at 30-35 °C. Reading at 24-48 until 72 h

According to ISO 11133 & USP Pharmacopeia

Microorganism

Burkholderia cepacea ATCC® 25416*Burkholderia cepacea* ATCC® 25608*Ps. aeruginosa* ATCC® 27853*Burkholderia cenocepacia* ATCC® BAA-245*Burkholderia multivorans* ATCC® BAA-247

Sterility Control

Add 5ml of the sample to 100 ml of TSB and to 100 ml Thioglycollate.

Incubation 48 hours at 30-35 °C and 48 hours at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

Growth

Good (≥50%) - Greenish-brown colonies with yellow halo

Good (≥50%) - Greenish-brown colonies with yellow halo

Inhibited

Good (≥50%). White colonies

Good (≥50%). White colonies surrounded by red zone

Bibliography

- COENYE, T., P. VANDAMME, J.R.W. Govan and J.J. LiPuma. 2001. Taxonomy and identification of the *Burkholderia cepacia* Complex. *J. Clin. Microbiol.* 39:10:3427-3436.
- ERAM, S.M., Q.B. NEJAD, G.R. Khatami, N. NAFISSI. 2004. Detection of *Burkholderia cepacia* Complex in Patients with cystic fibrosis. *Tanaffos* 3:9:47-52.
- GILLIGAN, P.H., P.A. GAGE, L.M. BRADSHAW, D.V. SCHIDLOW and LV DeCicco. 1985. Isolation medium for the recovery of *Pseudomonas cepacia* from respiratory Secretions of Patients with cystic fibrosis. *J. Clin. Microbiol.* 22:5-8.
- HENRY, d.Ä., M.E. CAMPBELL, J.J. LiPuma and D.P. SPEERT. 1997. Identification of *Burkholderia cepacia* Isolates from Patients with cystic fibrosis and use of a simple new selective medium. *J. Clin. Microbiol.* 35:3:614-619
- HENRY, d.Ä., M.E. CAMPBELL, C. McGimpsey, A. CLARKE, L. Louden, J.L. BURNS, M.H. ROE, P. VANDAMME and D.P. SPEERT. 1999. Comparison of isolation media for recovery from *Burkholderia cepacia* Complex Respiratory Secretions of Patients with cystic fibrosis. *J. Clin. Microbiol.* 37:4:1004-1007
- ISLES, A., I. MACLUSKY, M. COREY, R. GOLD, C. Prober, P. FLEMING, and H. LEVISON. 1984. *Pseudomonas cepacia* in cystic fibrosis Infection: An Emerging problem. *J. Pediatr.* 104:206-210.
- MILLER, JJ, 1996. *A guide to Specimen Management in Clinical Microbiology.* ASM Press., Washington D.C.
- WELCH, DF, M.J. Muszynski, C.H. PAI, M.J. MARCON, M.M. Hribar, P.H. Gilligan, J.M. Matsen, P.G. Ahlin, B.C.
- HILMAN and SA CHARTRAND. 1987. Selective and differential medium for recovery from the *Pseudomonas cepacia* respiratory tracts of Patients with cystic fibrosis. *J. Clin. Microbiol.* 25:1730-1734.
- ISO Standard 11133:2014 Microbiology of food, animal feed and water. Preparation, production, storage, and performance testing of culture media.
- USP 42 - NF 37 IS (2019) Chapter <60> Microbiological examination of nonsterile products: Test for *Burkholderia cepacia* complex. USP Corp. Inc. Rockville. MD. USA.