

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella Vi
Catalogue no.	TR1316
Package size	1 ml
Lot no.	2030323
Expiry date	2028-02
Intended use	Anti-Salmonella Vi records all <i>Salmonella</i> strains, as long as they contain antigen Vi in slide agglutination
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**