

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Coli III
Catalogue no.	TR2131-01
Package size	5 ml
Lot no.	6630823
Expiry date	2028-07
Intended use	Anti-Coli III records all <i>E. coli</i> strains, as long as they contain the antigens O 25:(K 11), O 78:(K 80), O 103:(K-), O 118:(K-), O 124:(K 72), O 145:(K-), O 157:(K-) and O 164:(K-).
Specificity	No agglutination with <i>E. coli</i> strains which are recorded by Anti-Coli I and Anti-Coli II.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**