

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh


Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:w</b>
Catalogue no.	TS1421
Package size	1 ml
<b>Lot no.</b>	<b>9271121</b>
<b>Expiry date</b>	<b>2026-10</b>
Intended use	Anti-Salmonella H:w records all <i>Salmonella</i> strains, as long as they contain the antigen combination H:ℓ,w in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
place, date

17.01.2022

  
signature Head of Quality Control