

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Coli I</b>
Catalogue no.	TS2111-01
Package size	5 ml
<b>Lot no.</b>	<b>4210523</b>
<b>Expiry date</b>	<b>2028-04</b>
Intended use	Anti-Coli I records all <i>E. coli</i> strains, as long as they contain the antigens O 26:(K 60), O 44:(K74), O 114:(K 90), O 125:(K 70), O 142:(K 86) and O 158:(K-).
Specificity	No agglutination with <i>E. coli</i> strains which are recorded by Anti-Coli II and Anti-Coli III.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2 °C...8 °C

**Dieses Dokument wurde elektronisch erstellt und ist ohne Unterschrift gültig.  
Die Chargenfreigabe erfolgte gemäß SOP QMV 7.**