

CERTIFICATE OF ANALYSIS

Product Name: Streptococcal Grouping - Extraction Reagent 3 Lot Number: C22299C

Product Code: PL.039 Expiration Date: 2026-01-31

Additional Information	
Storage Conditions:	

Instructions for Use: Refer to pro-lab.com/resources for current version

2-8°C

Quality Control Testing

Organism	Reaction*	Result
Streptococcus pyogenes ATCC 19615 (Group A)	Positive	PASS
Streptococcus agalactiae ATCC 12386 (Group B)	Positive	PASS
Streptococcus equisimilis ATCC 12388 (Group C)	Positive	PASS
Enterococcus faecalis ATCC 49477 (Group D)	Positive	PASS
Streptococcus sp. ATCC 12392 (Group F)	Positive	PASS
Streptococcus sp. ATCC 12394 (Group G)	Positive	PASS

*When used in conjunction with Extraction Reagents 1 and 2 to extract group specific carbohydrate antigens and tested against homologous latex reagents.

Release		
2024-03-12	Vampang	
Date	Quality Control	