

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 06-602-6300 | Revision: 10 |
| Remarks | | single use only |

Specification

| | |
|--|--|
| Product | Filter tip PP, premium surface, 0,1-10µl, short, 96 pcs/EcoRack, transparent,, mou. rings, np premium bio, sterile R |
| Dimensions | length: 31,30 mm (+/- 0,2 mm) inner diameter: 0,50 mm (+/- 0,2 mm) |
| Colour | transparent, |
| Graduation | mou. rings |
| Nominal Capacity | 0,1-10µl |
| Material Characteristics | polypropylene (PP), premium surface filter: high density polyethylene (PE) |
| Product Characteristics | biologically inert physiologically harmless no heavy metals according to EC/94/62 sterile R SAL 10-6 |
| Purity Grade | free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens and endotoxins Sterile R |
| Sterilisation Process | radiation sterilization SAL 10-6 |
| Shelf - Life | 60 months |
| Packaging | pcs per inner cs: 960 pcs per carton: 5.760 |
| Storage | dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable |
| Approved by nerbe plus quality department | |