

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:1
Catalogue no.	TR1437
Package size	1 ml
Lot no.	2540424
Expiry date	2029-03
Intended use	Anti-Salmonella H:1 records all <i>Salmonella</i> strains, as long as they contain antigen H:1,2, H:1,5, H:1,6, H:1,7. in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**