


Sabouraud Dextrose Broth EP/USP

N° Cat: 1205

N° Batch: 207282

Re-test Date: 7/2025

Regulations: USP / European Pharmacopoeia



Formula in g/l

Dextrose	20	Mixture of peptic digest of animal tissue and pancreatic digest of casein (1:1)	10
----------	----	---	----

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Batch pH
Sin restos	Polvo fino	Beige	Ámbar	5,6 ± 0,2	5,7

Microbiological test

According to European Pharmacopoeia, *Candida albicans*:
Incubation conditions: (30-35 °C / 3-5 days)

Rest of strains:
Incubation conditions: (30±2 °C / 18-48 h)

Microorganisms	Specification
<i>Candida albicans</i> ATCC 10231	Good growth, turbidity
<i>Aspergillus brasiliensis</i> ATCC 16404	Good growth, turbidity
<i>Candida albicans</i> ATCC 2091	Good growth, turbidity
<i>Escherichia coli</i> ATCC 25922	Partially inhibited
<i>Escherichia coli</i> ATCC 8739	Partially inhibited
<i>Lactobacillus rhamnosus</i> ATCC 9595	Good growth, turbidez
<i>Saccharomyces cerevisiae</i> ATCC 9763	Good growth, turbidity

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.