


GN Enrichment Broth (Hajna)

N° Cat: 1248

N° Batch: 204064

Re-test Date: 4/2025



Formula in g/l

Dextrose	1	D-mannitol	2
Potassium dihydrogen phosphate	1,5	Sodium chloride	5
Sodium citrate	5	Sodium deoxycholate	0,5
Tryptose	20	Potassium hydrogen phosphate	4

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Batch pH
w/o rest	Fine powder	Cream	Amber	7,0±0,2	7,0

Microbiological test

Incubation conditions: (35±2 °C / 6-24 h).

Microorganisms	Specification
Enterococcus faecalis ATCC 11700	Partially inhibited growth
Bacillus cereus ATCC 11778	Inhibited growth
Shigella flexneri ATCC 12022	Good growth
Salmonella typhimurium ATCC 14028	Good growth
Escherichia coli ATCC 25922	Good growth

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.