


Quality control certificate

Brain Heart Infusion Broth (BHI Broth)

N° Cat: 1400

N° Batch: 203034

Re-test Date: 3/2025



Formula in g/l

Dextrose	2	Disodium phosphate	2,5
Gelatin peptone	10	Sodium chloride	5
Heart infusion	10	Brain infusion	7,5

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Batch pH
w/o rests	Fine powder	Light toasted	Amber	7,4±0,2	7,5

Microbiological test

Incubation conditions: (35±2 °C / 18-24 h).

Microorganisms	Specification
Neisseria meningitidis ATCC 13090	Good growth with turbidity
Streptococcus pyogenes ATCC 19615	Good growth with turbidity
Streptococcus pneumoniae ATCC 6305	Good growth with turbidity

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.