

Quality control certificate

Haemophilus Test Medium

N° Cat: 1464

N° Batch: 208272

Re-test Date: 8/2025



Formula in g/l

Bacteriological agar	15	Hemin	0,015
Yeast extract	5	Mueller Hinton Base	21
NAD	0,015		

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Batch pH
w/o rests	Fine powder	Beige	Amber, slightly opalescent	7,3±0,1	7,4

Microbiological test

Incubation conditions: (35±2 °C / 16-18 h).

Microorganisms	Specification
Haemophilus influenzae ATCC 10211	Good growth

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.