

# Quality control certificate



## Haemophilus Test Medium

N° Cat: 1464

N° Lote: 105143

Re-test Date: 5/2024



### Formula in g/l

Bacteriological agar	15	Hemin	0,015
Yeast extract	5	Mueller Hinton Base	21
NAD	0,015		

### Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Lote pH
w/o rests	Fine powder	Beige	Amber, slightly opalescent	7,3±0,1	7,4

### Microbiological test

Incubation conditions: (35±2 °C / 16-18 h).

Microorganisms	Specification
Haemophilus influenzae ATCC 49247	Good growth

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.