

# Quality control certificate



## Kovac's Reagent

N° Cat: 5205

N° Lote: 820211

Re-test Date: 2/2022

### Formula in g/l

4-Dimethylaminobenzaldehyde (g)	50	Hydrochloric acid 37% (ml)	250
1-Butanol (ml)	750		

### Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Lote pH
w/o rests	Liquid	N/A	N/A	N/A	0,0

### Microbiological test

Incubation conditions: (35 °C / 24-48 h / anaerobic atmosphere).

Microorganisms	Characteristic reaction
Salmonella typhimurium ATCC 14028	Indol (-)
Escherichia coli ATCC 25922	Indole (+), Pink-red color

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.