

CERTIFICATE OF ANALYSIS

Product Name: Streptococcal Grouping Reagent Polyvalent Positive Control	Lot Number: C17533B
Product Code: PL.040	Expiration Date: 2021-02-28


Additional Information

Storage Conditions:	2-8°C
Precautions:	Refer to Instructions for Use

Quality Control Testing

Grouping Reagent	Reaction	Result
Streptococcal Grouping Reagent Latex Suspension Group A (PL.031)	Positive	PASS
Streptococcal Grouping Reagent Latex Suspension Group B (PL.032)	Positive	PASS
Streptococcal Grouping Reagent Latex Suspension Group C (PL.033)	Positive	PASS
Streptococcal Grouping Reagent Latex Suspension Group D (PL.034)	Positive	PASS
Streptococcal Grouping Reagent Latex Suspension Group F (PL.035)	Positive	PASS
Streptococcal Grouping Reagent Latex Suspension Group G (PL.036)	Positive	PASS

Release

2019 09 24	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill, Ontario,
 L4B 1K3
 support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
 21 Cypress Blvd., Suite 1070, Round Rock,
 Texas 78665-1034
 ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough, Wirral,
 Merseyside CH62 3QL
 uksupport@pro-lab.com