

### CERTIFICATE OF ANALYSIS

<b>Product Name:</b> Staph Negative Control Latex Reagent	<b>Lot Number:</b> C13078G	
<b>Product Code:</b> PL.086B	<b>Expiration Date:</b> 2017 05	


#### Additional Information

Storage Conditions:	2-8°C
Precautions:	Refer to Instructions for Use

#### Quality Control Testing

Appearance	Reaction*	Result
Homogenous suspension with no agglutination	N/A	Satisfactory
<b>Organisms</b>		
<i>Staphylococcus aureus</i> ATCC 25923	Negative	PASS
MRSA (Methicillin-resistant <i>Staphylococcus aureus</i> ) - Clinical Isolates	Negative	PASS
<i>Staphylococcus aureus</i> - Clinical Isolates	Negative	PASS
<i>Staphylococcus epidermidis</i> ATCC 12228	Negative	PASS
<i>Staphylococcus epidermidis</i> - Clinical Isolate	Negative	PASS
<i>Staphylococcus warneri</i> - Clinical Isolate	Negative	PASS
*Reaction within 20 seconds.		

#### Release

2015 10 09	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206  
20 Mural Street, Unit #4, Richmond Hill, Ontario,  
L4B 1K3  
support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450  
21 Cypress Blvd., Suite 1070, Round Rock,  
Texas 78665-1034  
ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614  
3 Bassendale Road, Bromborough, Wirral,  
Merseyside CH62 3QL  
uksupport@pro-lab.com