

CERTIFICATE OF ANALYSIS

Product Name: Staph Negative Control Latex Reagent	Lot Number: C16481F	
Product Code: PL.086B	Expiration Date: 2020-05-31	


Additional Information

Storage Conditions:	2-8°C
Precautions:	Refer to Instructions for Use

Quality Control Testing

Appearance	Reaction*	Result
Homogenous suspension with no agglutination	N/A	Satisfactory
Organisms		
<i>Staphylococcus aureus</i> ATCC 25923	Negative	PASS
MRSA (Methicillin-resistant <i>Staphylococcus aureus</i>) - Clinical Isolates	Negative	PASS
<i>Staphylococcus aureus</i> - Clinical Isolates	Negative	PASS
<i>Staphylococcus epidermidis</i> ATCC 12228	Negative	PASS
<i>Staphylococcus epidermidis</i> - Clinical Isolate	Negative	PASS
<i>Staphylococcus warneri</i> - Clinical Isolate	Negative	PASS
*Reaction within 20 seconds.		

Release

2018 10 04	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill, Ontario,
 L4B 1K3
 support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
 21 Cypress Blvd., Suite 1070, Round Rock,
 Texas 78665-1034
 ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough, Wirral,
 Merseyside CH62 3QL
 uksupport@pro-lab.com