

## **CERTIFICATE OF ANALYSIS**

Product Name:KOVACS REAGENTLot Number:K26934Product Code:PL.375Expiration Date:2019 11 30

Additional Information

Description: Pale green Concentrate Lot number: n/a Storage conditions: 15-25°C solution with characteristic Production date: 2018/11/30 Precautions: refer to SDS

odour.

## **Quality Control Testing**

Challenged with test organisms: Escherichia coli NCTC 10418

Enterobacter aerogenes NCTC 10006

## Test results:

Conforms to performance specifications when used in standard methods for detection of Indole production by coliform organisms in clinical, water and sewage samples.

## Comments:

We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

Release	
2018/11/30	matter.
Date	Quality Control

Canada № 800 268 2341 Fax 905 731 0206 20 Mural Street, Unit #4, Richmond Hill, Ontario, L4B 1K3 support@pro-lab.com U.S.A. 800 522 7740 Fax 800 332 0450 21 Cypress Blvd., Suite 1070, Round Rock, Texas 78665-1034 ussupport@pro-lab.com U.K. ≅ 0151 353 1613 Fax 0151 353 1614 3 Bassendale Road, Bromborough, Wirral, Merseyside CH62 3QL uksupport@pro-lab.com